

Stille  
Micro  
Instruments




**STILLE**

## Stille Micro Scissors

\* All patterns are supplied with round handles

Profile

Not shown at actual size



Art. no	Profile no.	Description
117-A110-15	1	Scissors with round handles, Sharp tips, 13 mm blades, Straight, 15 cm 6"
117-A110-18	1	Scissors with round handles, Sharp tips, 13 mm Blades, Straight, 18 cm 7"
117-A110-21	1	Scissors with round handles, Sharp tips, 14 mm Blades, Straight, 21 cm 8.25"
117-A111-15	2	Scissors with round handles, Blunt tips, 13 mm Blades, Curved, 15 cm 6"
117-A111-18	2	Scissors with round handles, Blunt tips, 13 mm Blades, Curved, 18 cm 7"
117-A111-21	2	Scissors with round handles, Blunt tips, 14 mm Blades, Curved, 21 cm 8.25"
117-A120-25	3	Scissors with round handles, Angled 25 degr, 11 mm blades, 18 cm 7"
117-A120-45	4	Scissors with round handles, Angled 45 degr, 11 mm blades, 18 cm 7"
117-A120-60	5	Scissors with round handles, Angled 60 degr, 11 mm blades, 18 cm 7"
117-A120-90	6	Scissors with round handles, Angled 90 degr, 11 mm blades, 18 cm 7"
117-A120-125	7	Scissors with round handles, Angled 125 degr, 11 mm blades, 18 cm 7"

## Stille Micro Needle holders

\* All patterns are supplied with round handles

Profile

Not shown at actual size



Art. no	Profile no.	Description
117-C110-15	1	Needle holders with round handles without lock, Diamond Jaws, Straight, 15 cm 6"
117-C110-18	1	Needle holders with round handles without lock, Diamond Jaws, Straight, 18 cm 7"
117-C110-21	1	Needle holders with round handles without lock, Diamond Jaws, Straight, 21 cm 8.25"
117-C111-15	2	Needle holders with round handles without lock, Diamond Jaws, Curved, 15 cm 6"
117-C111-18	2	Needle holders with round handles without lock, Diamond Jaws, Curved, 18 cm 7"
117-C111-21	2	Needle holders with round handles without lock, Diamond Jaws, Curved, 21 cm 8.25"
117-C120-15	1	Needle holders with round handles with lock, Diamond Jaws, Straight, 15 cm 6"
117-C120-18	1	Needle holders with round handles with lock, Diamond Jaws, Straight, 18 cm 7"
117-C120-21	1	Needle holders with round handles with lock, Diamond Jaws, Straight, 21 cm 8.25"
117-C121-15	2	Needle holders with round handles with lock, Diamond Jaws, Curved, 15 cm 6"
117-C121-18	2	Needle holders with round handles with lock, Diamond Jaws, Curved, 18 cm 7"
117-C121-21	2	Needle holders with round handles with lock, Diamond Jaws, Curved, 21 cm 8.25"

# Stille Micro Forceps


\* All patterns are supplied with round handles

## Profile

Not shown at actual size



Art. no	Profile no.	Description
117-B110-15	1	Forceps with tying platform, Straight, 15 cm 6"
117-B110-18	1	Forceps with tying platform, Straight, 18 cm 7"
117-B110-21	1	Forceps with tying platform, Straight, 21 cm 8.25"
117-B111-15	2	Forceps with tying platform, Curved, 15 cm 6"
117-B111-18	2	Forceps with tying platform, Curved, 18 cm 7"
117-B111-21	2	Forceps with tying platform, Curved, 21 cm 8.25"
117-B120-15	1	Forceps with tying platform, Diamond Jaws, Straight, 15 cm 6"
117-B120-18	1	Forceps with tying platform, Diamond Jaws, Straight, 18 cm 7"
117-B120-21	1	Forceps with tying platform, Diamond Jaws, Straight, 21 cm 8.25"
117-B121-15	2	Forceps with tying platform, Diamond Jaws, Curved, 15 cm 6"
117-B121-18	2	Forceps with tying platform, Diamond Jaws, Curved, 18 cm 7"
117-B121-21	2	Forceps with tying platform, Diamond Jaws, Curved, 21 cm 8.25"
117-B130-18	1	Forceps with tying platform, Diamond Jaws, Counter-balanced, Straight, 18 cm 7"
117-B130-21	1	Forceps with tying platform, Diamond Jaws, Counter-balanced, Straight, 21 cm 8.25"
117-B131-18	2	Forceps with tying platform, Diamond Jaws, Counter-balanced, Curved, 18 cm 7"
117-B131-21	2	Forceps with tying platform, Diamond Jaws, Counter-balanced, Curved, 21 cm 8.25"
117-B140-18	4	Forceps with ring tips 1.0x0.5 mm, Diamond Jaws, Counter-balanced, Straight, 18 cm 7"
117-B140-21	4	Forceps with ring tips 1.0x0.5 mm, Diamond Jaws, Counter-balanced, Straight, 21 cm 8.25"
117-B145-18	3	Forceps with ring tips 0.8x0.4 mm, Diamond Jaws, Counter-balanced, Straight, 18 cm 7"
117-B145-21	3	Forceps with ring tips 0.8x0.4 mm, Diamond Jaws, Counter balanced, Straight, 21 cm 8.25"



Serving surgery since 1841

## Stille Micro Instruments

Stille Micro Instruments are handmade instruments of the highest possible quality. The product line is focused mainly on plastic and vascular surgery. The black coating indicates that it is an original Stille Instrument.

All models are delivered with Stille's 10 year-warranty covering errors in material and workmanship.

Special models can be made on request.

### Stille Micro Scissors

The Stille Micro Scissors program covers straight and curved models in different lengths as well as vascular scissors with different blade angles. All models are made with the original Stille Super Cut blades for increased cutting performance in soft tissue.

The joint screw is designed to be removeable at service, and thereby recreate the same unique performance after service as the new instrument was delivered with.

### Stille Micro Forceps


All Stille Micro Forceps are designed to give the same unique feeling as the standard Stille forceps - a perfect spring tension in combination with perfect tip alignment and detail work. They are available with stainless steel tips or with Stille Diamond Jaws for increased durability.

The model program covers forceps with or without tying platform, ring tips and counter-balance, straight and curved models in different lengths.

### Stille Micro Needle holders

The Stille Micro Needle holders are designed to hold the finest sutures in a tight and secure grip with only a minimum of pressure from the surgeon's fingers.

All models are made with Stille Diamond Jaws. They can be delivered as straight or curved, with or without lock, in 3 different lengths. For suture sizes 6/0 and finer.



Stille Surgical AB

Hejargatan 8  
SE-632 29 Eskilstuna  
Sweden

Phone: +46--16-15 31 00  
Fax: +46-16-15 31 88  
e-mail: [info@stille.se](mailto:info@stille.se)  
[www.stille.se](http://www.stille.se)

Distributed by:

**TapMed** 

TapMed Medizintechnik Handel GmbH  
Gutshof 15 -17 • 34270 Schauenburg-Hoof

Telefon: 05601 - 9299-0  
Telefax: 05601 - 9299-44  
E-Mail: [info@tapmed.de](mailto:info@tapmed.de)  
Web site: [www.tapmed.de](http://www.tapmed.de)