



Titanium microsurgical instruments for vascular and nerve surgery



Excellent handling with premium-quality Biover Titanium Instruments


- + Non-magnetic titanium allows for precise handling of very fine sutures
- + Excellent, exact handling during surgery
- + Lightweight titanium reduces fatigue during long surgical procedures
- + Titanium has a non-reflective surface
- + Manufactured in cooperation with renowned surgeons
- + Made in Germany



1/1 (150 mm)

Titanium micro forceps for VASCULAR SURGERY

Ref. no.	Ø grip (mm)	Total length (mm)	tip (mm)	Platform
8038	8 mm	150 mm	0,2 mm straight	with tying platform
8039	8 mm	180 mm	0,2 mm straight	with tying platform



1/1 (150 mm)

Titanium micro needle holder for suture size 10/0, 11/0 and 12/0

Ref. no.	Ø grip (mm)	Total length (mm)	tip
8012	8 mm	150 mm	curved
8013	8 mm	180 mm	curved

Manufactured by:

Biover

Biover® AG

6052 Hergiswil | Switzerland

Tel. +41 (0)41 630 00 22 | biover@biover.ch

Distributed by:

VERKAUF

TapMed ✓

Medizintechnik Handels GmbH
Gewerbepark 10
34317 Habichtswald-Ehlen
GERMANY
Tel. +49 (0)5606 53060-0
Fax +49 (0)5606 53060-44

www.tapmed.de