


NEW

Titanium supermicro surgical instruments for supermicrosurgery




New: premium-quality Biover Supermicro Titanium Instruments

- + Supermicro instruments with handmade, precise, super fine tips for sutures up to 12/0
- + Non-magnetic titanium. The suture needle will not stick to the tips of the instruments.
- + Excellent, precise handling during supermicrosurgery
- + Lightweight titanium reduces fatigue during long surgical procedures
- + Titanium has a non-reflective surface
- + Manufactured in cooperation with renowned surgeons
- + Handmade in Germany



Titanium supermicro forceps for supermicrosurgery

Ref. no.	Ø grip (mm)	Total length (mm)	tip (mm)	Platform
8038.1	8 mm	150 mm	0,1 mm straight	with tying platform



Titanium supermicro needle holder for supermicrosurgery/suture size up to 12/0

Ref. no.	Ø grip (mm)	Total length (mm)	tip
8012.2	8 mm	150 mm	0.2 mm, curved

Manufactured by:

Biover

Biover® AG

6052 Hergiswil | Switzerland

Tel. +41 (0)41 630 00 22 | biover@biover.ch

Distributed by:

VERKAUF	TapMed <input checked="" type="checkbox"/>
	Medizintechnik Handels GmbH
	Gewerbepark 10
	34317 Habichtswald-Ehlen
	GERMANY
	Tel. +49 (0)5606 53060-0
	Fax +49 (0)5606 53060-44
	www.tapmed.de