

MICROSUTURE



MICROSUTURE is a nylon monofilament nonabsorbable, sterile surgical suture that is indicated for use in general soft tissue approximation and/or ligation, including use in microsurgical, cardiovascular, ophthalmic, and neurological procedures.

Features & Benefits

Characteristics

- Uniform diameter of the suture
- Maximum flexibility of the strand
- Consistent knot security
- Minimal memory



Ordering Information

Code	Size	Needle Type	Chord	Diameter	Curvature	Needle Length (arch)	Suture Length	Picture
GEM189BK	8-0 Nylon	taper	5 mm Chord	140 dia	3/8 curve	5.1 mm	15 cm/6" suture lgth	
GEM192BK	8-0 Nylon	taper	6 mm Chord	140 dia	3/8 curve	6.4 mm	15 cm/6" suture lgth	
GEM395BK	8-0 Nylon	tapersharp	6 mm Chord	140 dia	3/8 curve	6.4 mm	15 cm/6" suture lgth	
GEM177BK	9-0 Nylon	taper	3 mm Chord	100 dia	3/8 curve	3.8 mm	15 cm/6" suture lgth	
GEM180BK	9-0 Nylon	taper	4 mm Chord	100 dia	3/8 curve	5.1 mm	15 cm/6" suture lgth	
GEM188BK	9-0 Nylon	taper	5 mm Chord	140 dia	3/8 curve	5.1 mm	15 cm/6" suture lgth	
GEM386BK	9-0 Nylon	tapersharp	3 mm Chord	100 dia	3/8 curve	3.6 mm	15 cm/6" suture lgth	
GEM389BK	9-0 Nylon	tapersharp	5 mm Chord	100 dia	3/8 curve	5.1 mm	15 cm/6" suture lgth	
GEM171BK	10-0 Nylon	taper	3 mm Chord	70 dia	3/8 curve	3.8 mm	10 cm/4" suture lgth	
GEM385BK	10-0 Nylon	tapersharp	3 mm Chord	100 dia	3/8 curve	3.8 mm	10 cm/4" suture lgth	
GEM170BK	11-0 Nylon	taper	3 mm Chord	70 dia	3/8 curve	3.8 mm	10 cm/4" suture lgth	
GEM195BK	11-0 Nylon	taper	2 mm Chord	50 dia	3/8 curve	3.6 mm	10 cm/4" suture lgth	
GEM382BK	11-0 Nylon	tapersharp	3 mm Chord	70 dia	3/8 curve	3.8 mm	10 cm/4" suture lgth	

Black Needle

GEM192BKC	8-0 Nylon	taper	6 mm Chord	140 dia	3/8 curve	6.4 mm	15 cm/6" suture lgth	
GEM395BKC	8-0 Nylon	tapersharp	6 mm Chord	140 dia	3/8 curve	6.4 mm	15 cm/6" suture lgth	
GEM180BKC	9-0 Nylon	taper	5 mm Chord	100 dia	3/8 curve	5.1 mm	15 cm/6" suture lgth	
GEM389BKC	9-0 Nylon	tapersharp	5 mm Chord	100 dia	3/8 curve	5.1 mm	15 cm/6" suture lgth	

INDICATIONS FOR USE:

MICROSUTURE is indicated for use in general soft tissue approximation and/or ligation, including use in microsurgical, cardiovascular, ophthalmic, and neurological procedures.

Rx Only. For safe and proper use of this device, refer to the Instructions for Use.

Manufactured by Riverpoint Medical and Distributed by Synovis Micro Companies Alliance, Inc.